

Effect Diluent Low

BUFFERS & STABILIZERS

Sample Diluent & RF/HAMA Blocker

Cat. No. 7070

Lot No.: xxxxxxXxx

Date of manufacture: mm-yyyy

Expiry date: mm-yyyy (When stored at 2-8 °C in an unopened bottle)

Analysis: The diluent is inspected visually for appearance and the pH is measured.

Specifications: Appearance: Solution with slightly yellow to no distinct color
pH: 7.00 – 7.40

Results: Appearance: Passed
pH: Passed

Storage: Short-term storage: Ambient.
Long-term storage: 2-8 °C.

Performed by: **Date:** dd-MMM-yyyy **Signature:** XX

Approved by: **Date:** dd-MMM-yyyy **Signature:** XX

Our immunoassay solutions are environmentally conscious, creating a safer work environment while helping our customers fulfill significant regulatory requirements.
We call it ECO-TEK!



04-5-05-420 Ver. 4.1